

WHATMAN PAPER CHROMEDIA

Four new specialty
chromatographic papers
from Whatman.

GRADE

AC81

Acetylated cellulose
paper—12%

GRADE

AC82

Acetylated cellulose
paper—25%

For chromatography of
substances with very low water
solubility. Grade AC81 has an
acetyl content of approximately
12% and Grade AC82 an acetyl
content of approximately 25%.
Both grades are smooth surfaced
and have an excellent wet strength.

GRADE

SG81

Silica gel
loaded paper

A pure Whatman cellulose grade
containing a highly purified
silica gel adsorbent. Permits the
use of non polar solvents
and the separation of water
insoluble substances.

GRADE

AH81

Aluminum hydroxide
loaded paper

A pure Whatman cellulose paper
containing aluminum hydroxide.
Designed for the chromatographic
separation of substances
with very low water solubility.

For technical literature
request Bulletin PC-1.

Scientifica
P. O. BOX 1084, CLIFTON, NEW JERSEY

18. **Industrial Pharmacy** Section, American Pharmaceutical Assoc., 3rd annual eastern regional meeting, New York, N.Y. (H. Lieberman, Warner-Lambert Pharmaceutical Co., 170 Tabor Rd., Morris Plains, N.J.)

18-22. **Plastics**, 10th natl. exposition, Chicago, Ill. (J. Paluszek, G. M. Basford Co., 60 E. 42 St., New York 17)

18-22. **Radiation Accidents**, seminar on protection of the public, Geneva, Switzerland. (Office of the Director General, Food and Agriculture Organization, Viale delle Terme di Caracalla, Rome, Italy)

18-24. **Pan Indian Ocean Science** Assoc., New Delhi, India. (PIOSA, 39 Garden Rd., Karachi, Pakistan)

18-27. **Deans of Medical Schools**, conf., World Health Organization, Manila, Philippines. (WHO, Regional Committee for the Western Pacific, P.O. Box 2932, Manila)

19-21. **Building Research** Inst., fall conf., Washington, D.C. (BRI, 1725 De Sales St., NW, Washington)

19-21. **Contribution of Food Science** to Military Needs, Natick, Mass. (U.S. Army Natick Laboratories, Natick)

19-21. **American Physical Soc.**, Norman, Okla. (K. K. Darrow, 538 W. 120 St., New York 27)

19-21. **Stratosphere-Mesosphere** Structure conf., El Paso, Tex. (W. L. Webb, Schellenger Research Laboratories, Texas Western College, El Paso)

20-21. **Forest Biology**, 2nd conf., Mobile, Ala. (TAPPI, 360 Lexington Ave., New York 17)

21-22. **Pulmonary Diseases**, symp., Chapel Hill, N.C. (Div. of Pulmonary Diseases, Univ. of North Carolina School of Medicine, Chapel Hill)

21-23. **Forensic Physicians**, 10th conf., Prague, Czechoslovakia. (K. S. Lékarství, Tvrdeho ul 2a, Brno, Czechoslovakia)

21-23. **Oral Roentgenology**, research conf., Lincoln, Neb. (C. E. Crandell, Univ. of North Carolina School of Dentistry, Chapel Hill)

22-25. **Scientific Unions**, intern. council, Vienna, Austria. (ICSU 2, via Sebenico, Rome, Italy)

23. **American Translators Assoc.**, New York, N.Y. (D. S. Cunningham, German Dept., Rutgers Univ., 406 Penn St., Camden 2, N.J.)

24-27. **American Acad. for Cerebral Palsy**, Dallas, Tex. (J. D. Russ, 1520 Louisiana Ave., New Orleans, La.)

25-27. **Geological Soc. of America**, 76th meeting, New York, N.Y. (F. Betz, Jr., 419 W. 117 St., New York 27)


25-27. **Nuclear Electronics**, intern. symp., Paris, France. (French Society of Electronics Technicians, 10 avenue Pierre Larousse, Malakoff, Seine, France)

26-30. **Endocrinology**, 16th meeting, Paris, France. (Secrétariat du Service du Dr. Albeaux-Fernet, Hôpital Laënnec, 42, rue de Sévres, Paris 17)

29-30. **Biomagnetics**, intern. symp., Chicago, Ill. (J. F. Barnothy, Biomagnetic Research Foundation, 833 Lincoln St., Evanston, Ill.)

29-30. **American Mathematical Soc.**, Cleveland, Ohio. (AMS, 190 Hope St., Providence 6, R.I.)

31-1. **American College of Chest Physicians**, Portland, Ore. (M. Kornfeld, 112 E. Chestnut, Chicago 11, Ill.)



Yes, you can see what's going on inside the new **NALGENE® SEPARATORY FUNNEL . . .** right down to the stopcock. See the demarcation line clearly. Nalge's continuing research has produced the first polypropylene separatory funnel . . . unbreakable . . . translucent. Newest techniques permit blow molding the funnel in one continuous piece . . . nothing to impede the flow of the material. This separatory funnel (500 ml.) resists virtually all chemicals, and you can autoclave it, too. The non-sticking TEFLON* TFE plug provides a perfect seal without lubrication. No lubrication . . . no contamination. A polyethylene stopper won't freeze under any provocation. Every time a Nalgene separatory funnel is dropped, you're saving replacement dollars. It doesn't break! For more information write for Brochure M-563, Dept. 2133, The Nalge Co., Inc., Rochester 2, New York.

*Du Pont registered trademark.

 **NALGENE
LABWARE**
Leader in quality plastic labware since 1949