



FROM CORE LAB TO YOUR LAB.

Introducing the NEW easyCyte™ 8HT



- COMPLETE - the world's first end-to-end bench top flow cytometry solution.
- OPTIMIZED - multiparametric kits and assays.
- INTUITIVE - software and analysis packages.
- POWERFUL - innovative Guava® technology.
- AFFORDABLE - fits your bench top and budget.

Millipore brings the power of flow cytometry to your bench top.

Imagine running your most complex assays right on your bench top. No core labs. No outside expertise necessary. With Millipore's integrated instruments, assay kits, software and service, you have the end-to-end solution that makes it easy.

ADVANCING LIFE SCIENCE TOGETHER™ Research. Development. Production.

Millipore is a registered trademark of Millipore Corporation.
Guava is a registered trademark of Guava Technologies, Inc.
Advancing Life Science Together and the Millipore logo are trademarks of Millipore Corporation.
easyCyte is a trademark of Guava Technologies, Inc.
©2009 Millipore Corporation. All rights reserved.



Make the move to Guava at
www.millipore.com/flowcytometry5

Submission
deadline
August 1

Rewarding brilliance since 1995.



The GE & Science Prize for Young Life Scientists. Because brilliant ideas build better realities.

Imagine standing on the podium at the Grand Hotel in Stockholm, making your acceptance speech. Imagine joining the ranks of those published in *Science* magazine and having your essay on your work in molecular biology read by your peers around the world. Imagine taking part in a seminar with the other Prize winners and Nobel Prize laureates and discussing your work with leaders in the field. Imagine what you could do with the \$25,000 prize money. Imagine what a brilliant start to your career, and where it could lead. Now stop imagining. If you were awarded your Ph.D in 2008, submit your 1000-word essay by August 1 and make it reality.

Want to build a better reality? Go to www.gelifesciences.com/science



GE & Science
Prize for Young
Life Scientists



imagination at work



* For the purpose of this prize, molecular biology is defined as "that part of biology which attempts to interpret biological events in terms of the physico-chemical properties of molecules in a cell".

(*McGraw-Hill Dictionary of Scientific and Technical Terms, 4th Edition*).

GE Healthcare Bio-Sciences AB, a General Electric Company.
Björkgatan 30, 751 84 Uppsala, Sweden.
© 2009 General Electric Company
- All rights reserved.

Biacore systems



from inspiration
...to publication

Highest quality, information-rich interaction data from Biacore™ systems deepen your understanding of molecular mechanisms and interaction pathways and enable you to add function to structure.

Select the perfect solution for your application and draw conclusions with confidence – from the company that continues to set the standard for label-free protein interaction analysis.

For further information or register to have one of our scientific experts contact you, please visit www.gelifesciences.com/biacore-science



Biacore T100
unmatched performance




Biacore X100
ready to run research system




Biacore Flexchip
array-based comparative profiling




imagination at work



Immunohistochemical analysis of paraffin-embedded human astrocytoma using **VEGF Receptor 2 (55B11) Rabbit mAb #2479**.



Immunohistochemical analysis of paraffin-embedded mouse small intestine using **MMP7 (D4H5) Rabbit mAb #3801**.




The Highest Quality Antibodies for the Study of

Angiogenesis

...from Cell Signaling Technology®

Unparalleled product quality, validation and technical support.



Confocal immunofluorescent analysis of HeLa cells, treated with either 10 μ M MG132 (left) or 10 μ M MG132 and 1 mM DMOG (right), using **Hydroxy-HIF-1 α (Pro564) (D43B5) Rabbit mAb #3434** (green). Actin filaments have been labeled using DY-554 phalloidin (red).

- ❑ Cell Signaling Technology's innovative products offer unsurpassed sensitivity, specificity and performance.
- ❑ Extensive in-house validation means that optimization is not left up to you, the user.
- ❑ Technical support provided by the same scientists who produce and validate the products translates into a thorough, fast and accurate response.

for quality products you can trust...

www.cellsignal.com




Cell Signaling

TECHNOLOGY®





Advanced Transfection Technologies!

- ▶ *High efficiency, low cytotoxicity*
- ▶ *Simple, serum-compatible protocol*



Xfect performs 40-fold better than a leading competitor (Product LX) on Jurkat cells, while Xfect Stem performs 6-fold better on mouse ES cells.

Sign up for a free sample.
www.clontech.com/xfect2



Effective transfection of embryonic stem cells with Xfect Stem. ES-E14TG2a mES cells were transfected with pDsRed-Express2 using Xfect Stem. After 48 hr, transfection efficiency was assessed via flow cytometry (A) as well as by phase and fluorescence microscopy (B & C). As quantified by flow cytometric analysis, 93.9% of the stem cells were positive for DsRed-Express2 expression.

Clontech Laboratories, Inc.
 A Takara Bio Company
www.clontech.com

United States/Canada: +1.800.662.2566 • Asia Pacific: +1.650.919.7300 • Europe: +33.(0)1.3904.6880 • Japan: +81.(0)77.543.6116
 For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Clontech, the Clontech logo, and all other trademarks are the property of Clontech unless noted otherwise. ©2009 Clontech Laboratories, Inc.

AD963177 US





FROM START TO FINISH, ENDNOTE X3 KEEPS YOUR RESEARCH ON THE FAST TRACK.

With its world-class speed and track record, EndNote X3 puts your publishing in the lead and keeps you there. Out of the blocks with new features like faster start up times and seasoned favorites like our Cite While You Write™ technology, EndNote X3 is moving fast to deliver the world's most complete bibliographic solution.

It doesn't stop there. EndNote X3 is also meters ahead with a personal EndNote Web account that not only transfers groups swiftly between desktop and Web but also organizes your own publication list for the ResearcherID author community. You'll find new speed in Cite While You Write for Apple® Pages '09, Microsoft® Word, and now OpenOffice.org Writer 3 for Windows. EndNote X3 even sports new options for chemistry styles, multiple bibliographies within a Word document, and more.

Keep the lead, right to the finish line. Put EndNote X3 on your reference team today!

800-722-1227 • 760-438-5526 • rs.info@thomson.com



Download your free demo or buy online today
www.endnote.com




Switch to Superior High Fidelity PCR

PrimeSTAR[®]

Top Five Reasons to Switch to PrimeSTAR[®]:

- ▶ **Higher Accuracy:** A strong exonuclease activity results in an extremely low error rate, with only 15 errors per 480,000 bp on a GC-rich template.
- ▶ **Higher Efficiency:** Higher than Taq Polymerase.
- ▶ **Robust Amplifications:** A single PCR cycling protocol can be used to amplify products of varying sizes.
- ▶ **Greater Specificity:** The antibody mediated Hot Start formulation prevents false initiation events during reaction assembly.
- ▶ **Excellent with GC-Rich Targets:** Robust performance with GC-rich templates using the GC buffer formulation.

Amplification Efficiency of a 2 kb Human Genomic DNA Fragment (DCLRE1A). Excellent sensitivity and yield are observed when amplifying with PrimeSTAR[®]. Quantities: Lane 1: 0 ng (dH₂O), Lane 2: 100 pg, Lane 3: 1 ng, Lane 4: 10 ng, Lane 5: 100 ng.



PrimeSTAR[®] is a registered trademark of Takara Bio Inc. Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights (such as 5' Nuclease Process patent rights) are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA. Takara Bio's Hot-Start PCR-Related products are licensed under U.S. Patent 5,338,671 and 5,587,287 and corresponding patents in other countries.

Takara

For more information
www.takara-bio.com

Japan:
 Takara Bio Inc.
 +81 77 543 7247
www.takara-bio.com

USA:
 Takara Bio USA
 A Division of Clontech Laboratories, Inc
 888-251-6618
www.takarabio.com

Europe:
 Takara Bio Europe S.A.S.
 +33 1 3904 6880
www.takara-bio.eu

China:
 Takara Biotechnology
 (Dalian) Co., Ltd.
 +86 411 8764 1681
www.takara.com.cn

Korea:
 Takara Korea
 Biomedical Inc.
 +82 2 2081 2525
www.takara.co.kr



Imaging Systems

The Leica MM AF imaging systems combine Leica's microscopy with MDS Analytical Technologies' MetaMorph software to address a range of applications in wide-field imaging in life sciences research. The tightly integrated systems are based on the ease-of-use and optical performance of the Leica microscope range, giving researchers an efficient platform for many types of imaging experiments. An integrated journaling capability allows flexible customization of the platform for virtually any imaging application.

Leica Microsystems

For information +49-(0)-6441-29-2550
www.leica-microsystems.com

Petri Dish Filler

A compact, high-volume media preparation system provides a fast, automated method for filling Petri dishes. The Integra Mediajet is designed for high-volume microbiological media users who handle up to 20,000 plates per week. All functions can be controlled via a full-size screen. Once the dish-filling process has started the system continues the cycle through to completion without the need for input or monitoring, allowing operators to work on other tasks. The system is designed for reliable and independent operation, as the dishes are mechanically guided through the process and monitored by a set of sensors during the filling process. Typical variations in the diameter or shape of disposable plastic dishes, as well as inconsistencies in the rims or lids, are also compensated for by the system, eliminating the potential for dish jams.

Scientific Laboratory Supplies

For information +44-(0)-1159-821-111
www.scientific-labs.com

Low Particulate Bottles

Nalgene Certified Clean PETG Containers are suitable for pharmaceutical and biotechnology applications requiring the processing and storage of bulk intermediates and critical reagents such as vaccine and protein therapeutic preparations. The containers provide excellent gas barrier properties to minimize pH shifts. The resins used to produce these containers are free of animal-derived components. The bottle with its closure system is guaranteed leakproof, making it suitable for use in transportation.

Thermo Fisher Scientific

For information 508-742-5254
www.thermofisher.com

Ion Analysis

Two new kits are available for anion and cation analysis on the P/ACE MDQ and ProteomeLab PA 800 Systems. Capillary electrophoresis analysis of ions performed using these kits delivers efficient, rapid separations with small sample requirements and short preparation times. The kits incorporate a unique dynamic coating technology that produces consistent migration times. In drug discovery applications, the kits provide robust results for counter-ion analysis, an important part of purity determination. Compound libraries can be processed for counter-ion content even if solubility data are not available. The kits allow analysis of compounds over a wide range of

polarities for applications with diverse sample matrices. These kits can also be used in the analysis of foods, beverages, and industrial products for detection and quantitation of organic acids, inorganic anions and cations, and aliphatic amines.

Beckman Coulter

For information 800-742-2345
www.beckmancoulter.com

Cap and Tubes Strips

Masterclear Cap Strips and real-time PCR Tube Strips optimize light transmission and bring the benefits of white reflective wells to the polymerase chain reaction (PCR) tube format. The white wells of the real-time PCR Tube Strips significantly increase the fluorescent signal compared with clear or frosted wells, and reduce the interference of the thermal block. The result is improved sensitivity and reproducibility in real-time PCR experiments, which may enable the number of replicates to be reduced. The inverted dome design of Masterclear Cap Strips prevents scratching or contamination of the optical surface, while the extremely thin transmission windows ensure transmission values greater than 90 percent for wavelengths above 500 nm, with low variance. The increase in fluorescence signal can thus be distinguished from background noise at an earlier time point. The products provide particular advantages in detection of low-copy genes, when using reagents with weak fluorescent signals, and in low-volume real-time PCR assays.

Eppendorf

For information +49-40-538-01-640
www.eppendorf.com/consumables

Adhesive Films

Pattern Adhesive Films provide high-integrity seals for 96-well and 384-well microplates. They have adhesive printed in a pattern that exactly matches the shape of a microplate well, so no adhesive comes into contact with samples or solvent. In contrast to traditional adhesive seals, piercing a Pattern Adhesive Film with a single or multichannel pipettor or robotic probe does not risk contaminating a sample with residual adhesive. Made from 50- μ m thick PTFE film, they are tough yet easily pierced. They offer excellent optical reading properties.

Porvair Sciences

For information +44-1372-824290
www.porvair-sciences.com

Electronically submit your new product description or product literature information! Go to www.sciencemag.org/products/newproducts.dtl for more information.

Newly offered instrumentation, apparatus, and laboratory materials of interest to researchers in all disciplines in academic, industrial, and governmental organizations are featured in this space. Emphasis is given to purpose, chief characteristics, and availability of products and materials. Endorsement by *Science* or AAAS of any products or materials mentioned is not implied. Additional information may be obtained from the manufacturer or supplier.

The Proof of Pluripotency.

Now You *Know* Your ES and iPS Cells Are Pluripotent.

Because it's all about pluripotency. Stem cell researchers know all too well that pluripotency is far from a reliable or steady state. Any number of factors or signaling cascades can cause cells to differentiate. If that happens, you've lost time, money and critical research data. That's why we've introduced the Stemgent Antibody Collection.

- ❑ **Accurate characterization is the cornerstone of results you can trust.**
- ❑ **Our antibodies are highly specific for the most universally accepted pluripotency markers.**
- ❑ **Validating on ES cells provides more significant results.**

MORE www.stemgent.com/ab2

High Performance Antibody Collection From Stemgent.

We've done the screening for you. Multiple antibodies from multiple sources are validated for specificity, sensitivity, and reproducibility to their corresponding protein targets. And Stemgent takes the time to verify the stem cell relevance of each antibody by performing ICC and FC analysis on embryonic stem cells.


- ❑ **The entire collection is Application Tested on stem cells.**
- ❑ **High performance antibodies to ensure success.**
- ❑ **Quickly becoming the most relevant stem cell antibody collection available.**




MORE www.stemgent.com/ab4

Our purpose is to move stem cell science forward.

The best way to do that is to make the research of world-class scientists available to you in ways that are reliable and useful. In that way, your work benefits from their work. And the future of stem cell research is accelerated.



A and B. ICC analysis of the Stemgent® Purified Rabbit anti-Mouse/Human Sox2 Antibody on mouse ES cell line R1 (image A), on human ES cell line H1 (image B)



C and D. ICC analysis of the Stemgent® Affinity Purified anti-Human TRA-1-60 on human ES cell line H1 (image C), ICC analysis of the Stemgent® Affinity Purified anti-Human SSEA-3 on human ES cell line H1 (image D)

Save Time and Resources by Using Pre-Screened Antibodies.

Validating antibodies from a choice of vendors is often a complex, time-consuming project. We are here to help. We invite you to be in direct contact with a fully qualified Stemgent antibody specialist who can help guide you through our unique collection, and choose the right antibodies for you.

- ❑ **Intracellular and cell surface markers**
- ❑ **Affinity purified and conjugated formats**
- ❑ **Isotype controls**

MORE www.stemgent.com/ab6

Come to www.stemgent.com/sci11 for more information. Or call 877-228-9783 (toll-free) or +1-617-245-0098 (international).



BOSTON | www.stemgent.com/sci11 | SAN DIEGO

MORE www.stemgent.com/xxx

About Stemgent Information Links. Get to the precise information you need faster. When you see one of these link bars, go to Stemgent.com and add the /xxxx extension. It's a direct link to the info you need.

NOTE: Stemgent acknowledges that the TRUE test of pluripotency is differentiation into all 3 germ layers. In many cases it is more important to get the most information you can get FAST. The Stemgent® Antibody Collection gets you the results you need on your schedule.

© 2009 by Stemgent, Inc. Stemgent is a registered trademark and Reprogramming the Reagent, and the What's Next logotype are trademarks of Stemgent, Inc.

Everyday needs. Everyday



answers™

 invitrogen™

Easily find your favorite brands of reagents, kits, and benchtop devices at www.invitrogen.com/everyday.

Ambion®
RNA and RNA sample prep

Dynal®
Magnetic beads for separation

Gateway®
Cloning and protein expression

GIBCO®
Cell culture

Lipofectamine™
Transfection reagents

Molecular Probes®
Fluorescent dyes and probes

Novex®
Protein separation and blotting

SuperScript®
cDNA synthesis reagents

TaqMan®
Real-time PCR

TOPO®
PCR cloning

For research use only. Not intended for any animal or human therapeutic or diagnostic use, unless otherwise stated.

© 2009 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. TaqMan is a registered trademark of Roche Molecular Systems, Inc. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. A-085988 0609