

biomapping



Bioextract.

Overcome protein sample complexity –
Separate with Sigma's Seppro[®]
Depletion Technology.

Sigma's Seppro Depletion Technology enables removal of interfering highly abundant proteins from a variety of biological samples using the affinity of avian polyclonal IgY antibodies. This removal un.masks low abundance proteins of interest for further investigation.

wherebiobegins.com/biomapping

SIGMA-ALDRICH[®]

SIGMA[®] *Where bio begins*[™]
Life Science

Seppro is a registered trademark of Sigma-Aldrich Biotechnology LP and Sigma-Aldrich Co.