

Now There's an Easier Path to Extracting RNA and Protein

CellAmp®

CellAmp® Direct RNA Prep Kit for RT-PCR (for Real Time) & Protein Analysis

CellAmp® Direct RNA Prep Kit for RT-PCR (Real Time) & Protein Analysis is designed for the preparation of cell lysates that can be used for either one or two step RT-qPCR mRNA expression analysis or Western blot protein expression analysis of cultured animal cells. By following a simple protocol (Figure 1), a lysate can be prepared in as little as 10 minutes. In addition, the prepared cell lysate

may be used as either a template for RT-qPCR or as a sample in a regular Western blot analysis without further nucleic acid purification or protein extraction.

- ✓ Quick and Simple Preparation of RNA
- ✓ Gene Expression Analysis Performed in 2 hours

Features

- **Reproducible, High Quality Data**
- **Lysates can be used for both mRNA and Protein Expression Analysis** using either RT-qPCR or Western Blotting
- **Prepare Cell Lysates in only 10 Minutes**

The use of these products and technologies may be governed by one or more licensing agreements. To view specific and current licensing information for any Takara Bio products, please consult the product webpage on our website at www.takara-bio.us.




Figure 1: CellAmp® Direct RNA Prep Kit Protocol

Takara

For more information
www.takara-bio.com

USA:
Clontech Laboratories
www.takara-bio.us
800-662-2566

Visit our Website Today!
www.takara-bio.com