

Rethink Western blotting. Take protein detection to new dimensions.

Biology isn't flat—neither are your Westerns.

Meet the fast, versatile SNAP i.d.® 2.0 system, to fully exploit three-dimensional reagent distribution. Unlike conventional Western blotting, where diffusion is the primary means of reagent transport, the SNAP i.d.® 2.0 system applies a vacuum to actively drive antibodies and buffers right through the membrane. This advanced technology promotes antigen binding and thorough washing, enabling you to better optimize your Western blotting conditions.

www.emdmillipore.com/SNAP



EMD Millipore is a division of Merck KGaA, Darmstadt, Germany