

Equipment for

**FLUORESCENCE
MICROSCOPY**

and

**ULTRAVIOLET
PHOTOMICROGRAPHY**

A NEW, high-intensity illuminating apparatus, applicable to any standard microscope, for examination of opaque or transparent specimens under fluorescent light. No perceptible loss of brilliancy with high magnifications.

Biological specimens may be examined without staining, or after a brief treatment with fluorochrome. Different animal and plant tissues exhibit characteristic colorations which make possible hitherto difficult differentiations. Adaptable for observation of living organs in situ, without interference with function or structure. Also recommended for paper, oils, rubber, pigments, textile fibers, cellulose, etc.

Same source of light may be used for photomicrography, as it gives an intense, steady light, covering the entire spectral range. Clearly defined lines at 3650 and 5460 Å, thus permitting full advantage to be taken of the increased resolution and differentiation obtainable with the ultraviolet method.

Particulars on request

Pfaltz & Bauer, Inc.
EMPIRE STATE BUILDING, NEW YORK

*The extreme accuracy
of precision machined*

HOLMES Projectors


assures rock steady brilliant
pictures at all times.

Accurate projector speed.

Clear toned and powerful
sound reproduction

Minimum annoyance to audience
because of film breakage, illumina-
tion failures or slow threading.

Built with the mechanical
precision of largest
projectors, yet costs no
more than portables.



16 mm HOLMES
Ball Bearing—All
shaft driven—All
sprocket.

Many exclusive features found only
in a HOLMES is the reason it is
favored by professional operators
for portable uses.

Write for latest catalog and prices.

**HOLMES 16 mm
35 mm PROJECTORS**


Choice of Professional Operators Everywhere

HOLMES PROJECTOR CO.

1821 Orchard St., Chicago.

ZEISS

COMPLETE EQUIPMENT FOR USE IN QUANTITATIVE SPECTRUM ANALYSIS



QUARTZ SPECTROGRAPH QU 24 possesses an optical system of high luminosity. Maximum effective aperture ratio of $F/11.5$. Exceptional definition combined with flat spectrum. Uniform resolution over whole ratio from 3,000 to 5,800 Au ; lines with a difference of 0.348 Au are sharply separated. Uses ordinary plates of standard thickness. Rigidly constructed, simple in operation, adjustments unnecessary.

FEUSSNER SPARKING APPARATUS (Patented). A rotating spark gap driven by a synchronizing motor assures perfect constancy of spark necessary for accurate quantitative analysis.

SPECTRUM LINE PHOTOMETER permitting accurate and reliable objective photometric evaluation of spectrograms for quantitative analysis.

SPECTRUM PROJECTOR supplies a sharply defined image of the spectrogram at a magnification suitable for visual evaluation.

CARL ZEISS, INC., 485 FIFTH AVE., NEW YORK
728 So. Hill Street, Los Angeles