

Class E Spectroscopic Plates

THIS RECENT ADDITION to the series of *Eastman Spectroscopic Plates* has its highest sensitivity through the region 6000Å. to 7000Å., with a maximum at 6500 Å. Class E sensitizing can be supplied with Types I, II, III, IV, V, and 144 emulsions.

Astronomers have found these plates valuable in the photography of faint red stars, and they have been successfully used by spectroscopists for photographing weak lines in the red region of the spectrum.

Additional information concerning photography in the H-alpha region will be supplied on request.

EASTMAN KODAK COMPANY
Research Laboratories **ROCHESTER, N. Y.**

Sulfonphthalein Indicators

Savings in manufacturing costs permit lower prices to you on GUARANTEED indicators!

Thymol Blue, per gram \$.45 in 10-gram vials.

Order direct or from your dealer

The Synthetical Laboratories

5556 Ardmore Ave.

Chicago, Ill.

The LaMotte Falling Drop Densimeter

for determining

Blood Proteins and Specific Gravity of Body Fluids

(after Barbour and Hamilton)

Blood proteins determined with one drop of plasma. Blood and plasma specific gravity changes can be followed minute by minute.

Method is accurate to 0.0001 specific gravity.

Requires less than a minute to determine specific gravity.

Write for new illustrated descriptive bulletin.

LaMotte Chemical Products Co.

Dept. "H"

Towson, Baltimore, Md.

Bargains in New and Used Glasses! All makes,



Size & powers of nature, sport & hunting glasses, \$6. upward. Telescopes, spotting & riflescopes. Microscopes, magnifiers & compasses \$1. upward. Repair work & goods guaranteed.

J. ALDEN LORING, OWEGO, N. Y.

Catalog free, Box B.

SPRAGUE - DAWLEY, INC.

Pioneers in development of the standard laboratory rat

Madison, Wisconsin

Illuminating Pocket Magnifier

CARL ZEISS
JENA



ZEISS

Repair Service

Our stock of CARL ZEISS optical instruments, accessories and repair parts on hand in New York is sufficient for a considerable time. Purchasers and owners of Zeiss equipment can be assured of uninterrupted service in our large repair shop at our main office in New York, or in the repair shop of our Los Angeles branch.

THE ZEISS Illuminating Pocket Magnifier makes possible accurate observation under poor lighting conditions by providing its own source of illumination. It is useful for examining the surfaces of scientific, mineralogical or geological specimens, as well as the interior of cavities, &c. Four interchangeable slip-on mounts for the aplanatic magnifying lenses provide a range of 3x, 6x, 8x or 10x magnification.

Leaflet MED 333 upon request.

CARL ZEISS, INC.

485 Fifth Avenue
NEW YORK

728 So. Hill Street
LOS ANGELES

CARL ZEISS
JENA