

Cobalt-57  
 Cobalt-58  
 Cobalt-60  
 Copper-64  
 Gadolinium-153  
 Gallium-67  
 Germanium-68  
 Gold-195  
 Gold-198  
 Hydrogen-3  
 (Tritium)  
 Indium-111  
 Indium-111  
 Oxine  
 Indium-114m  
 Iodine-125  
 Iodine-129  
 Iodine-131  
 Iron-55  
 Iron-59  
 Krypton-85  
 Manganese-54  
 Mercury-203  
 Molybdenum-99  
 Nickel-63  
 Phosphorus-32  
 Phosphorus-33  
 Potassium-42  
 Promethium-147  
 Rubidium-86  
 Ruthenium-103  
 Samarium-151  
 Scandium-46  
 Selenium-75

Silver-110m  
 Sodium-22  
 Sodium-24  
 Strontium-85  
 Strontium-90  
 Sulfur-35  
 Technetium-99  
 Thallium-204  
 Thulium-171  
 Tin-113  
 Tin-119m  
 Tin-121m  
 Titanium-44  
 Ytterbium-169  
 Zinc-65  
 (and two dozen  
 others available  
 periodically)

Antimony-125  
 Barium-133  
 Bismuth-207  
 Bromine-82  
 Cadmium-109  
 Cadmium-115m  
 Calcium-45  
 Carbon-14  
 Cerium-141  
 Cesium-134  
 Cesium-137  
 Chlorine-36  
 Chromium-51

**Nobody offers  
 more radionuclides  
 or faster response than  
 we do. Call 800-225-1572.  
 We'll have the nuclide you need,  
 packaged the way you want it,  
 on its way to you tomorrow.**

Not for use in humans or clinical diagnosis.

New England Nuclear  
 549 Albany Street, Boston, MA 02118  
 Call toll free: 800-225-1572, Telex: 94-0996  
 Mass and Internat'l: 617-482-9595  
 Europe: NEN Chemicals GmbH, D-6072, W. Germany  
 Postfach 401240, Tel. (06103) 803-0, Telex 4-17993 NEN D  
 NEN Canada: 2453 46th Avenue, Lachine, Que H8T 3C9  
 Tel. 514-636-4971, Telex 05-821808 © 1982 NEN


**New England Nuclear**  
 a Du Pont company


Circle No. 201 on Readers' Service Card

## TINY FLOW METERING PUMPS WITH AC MOTOR DRIVE

Those handy little RH pump heads for  
 piddling flow rates now come packaged  
 with their own synchronous AC motor  
 drive. Each has an RH miniature pump head and  
 totally enclosed three-speed motor drive  
 assembly. This provides a choice of 150, 300 or  
 600 piston strokes per minute by means of a safe  
 and simple belt change arrangement. Motor can  
 be stalled indefinitely without overheating!

Pumps plug into standard 110-volt 60-Hz three-prong  
 sockets. Are steplessly adjustable even while running  
 from zero to maximum flow rate by means of calibrated  
 control ring. Are reversible while operating. And have  
 the exceptional chemical resistance of alumina ceramic  
 and fluorocarbon in the wetted stream. Compression  
 fittings at the ports accept 1/4 OD standard laboratory  
 tubing. Pressures to 50 psi at low flows, to 20 psi at  
 any flow.

Repeat Accuracy to better than 1%!



Model RHSY0CKC \$380  
 Max flow rate: 8.5 ml/min at 150 SPM, 17 ml/min  
 at 300 SPM, 34 ml/min at 600 SPM  
 Model RHSY1CKC \$380  
 Max flow rate: 15 ml/min at 150 SPM, 30 ml/min at  
 300 SPM, 60 ml/min at 600 SPM  
 Available from stock for immediate delivery. Write or  
 call for Catalog RP401: the FMI LAB PUMP STORY.


**FLUID METERING, INC.**  
 29 ORCHARD STREET, OYSTER BAY, N.Y. 11771 • (516) 922-6050

Circle No. 168 on Readers' Service Card

**NEW! EXCLUSIVE!**

# $\alpha$ -EGF

**HPLC PURIFIED  
EPIDERMAL GROWTH FACTOR**

A HPLC purified homogenous form of Epidermal  
 Growth Factor.  $\alpha$ -EGF is many times more  
 potent than other Epidermal Growth Factor  
 preparations. Biomedical Technologies is a  
 primary manufacturer of Growth Factor  
 Products:

EGF (P-10, DEAE)	Hybridoma Cell Growth Supplement
Endothelial Mitogen	2.5S NGF
FGF (Pituitary)	7S NGF
FGF (Brain)	PdGF (>50,000units/mg)
Fibronectins	$\alpha$ -Thrombin

For Research Use Only

For more information, please call or write:

**bti**

Biomedical Technologies Inc.  
 22 Thomdike Street  
 Cambridge, MA 02141  
 (617) 492-6200

Circle No. 214 on Readers' Service Card