

What can **RabMAbs[®]** do for you?

#3	Excellent for IHC
#4	High Affinity
#5	High Specificity
#6	Novel Epitopes

RabMAbs[®] consistently show higher affinity, with K_D values which can often reach the picomolar level (10^{-12} M), effectively eliminating the need for further affinity maturation.

LOW AFFINITY 10^{-5} 10^{-12} 10^{-15} HIGH AFFINITY

Antibody affinity average value for a typical RabMAb from analysis of over 850 RabMAbs[®] (1).

1. Pope ME, *J Immunol Methods*. 2009 Feb 28;341(1-2):86-96.



Rabbit Monoclonal Antibodies (RabMAbs[®]) offer multiple advantages to bring you the highest quality antibody possible.

Discover more at abcam.com/RabMAbs

