



Less is more.

Less Work.
More Dependability.

They're ready-to-use kits with everything you need except your sample. Room temperature set-up, no salting out, no second incubation. We've pre-reacted the first antibody with the second antibody to save you time.

And we've given you NEN quality to assure you of dependable assays — accurate and reproducible. Sensitivity reaches the femtomole range with the acetylating reagents provided. A tritium-labeled marker is included to monitor recovery. The iodinated tracer is exceptionally stable.

These are some of the reasons to choose NEN's Cyclic AMP and GMP [¹²⁵I] RIA kits, the world's most widely used for their convenience and quality. Send for the whole story in our technical brochure.

Not for use in humans or clinical diagnosis.

NEN New England Nuclear

549 Albany Street, Boston, Mass. 02118
Call toll-free: 800-225-1572
(In Massachusetts and International: 617-482-9595)
NEN Chemicals GmbH: D-6072 Dreieich, W. Germany
Postfach 401240, Telephone (06103) 85034, Telex 4-17993 NEN D
NEN Canada Ltd., 2453 46th Avenue, Lachine, Que H8T 3C9.
Telephone 514-636-4971, Telex 05-821808

The New Drummond Digital Microdispenser.

You can give better than lip service to volumetric accuracy.

Get accuracy and reproducibility of better than $\pm 1\%$. With the all-new Digital Microdispenser from Drummond® — the originators of positive displacement pipetting. The unique Drummond dual-check on accuracy assures that the liquid in the capillary equals the reading on the handle. And we've made it as easy as possible to get repeated precision delivery (with one hand). And as hard as possible to make a mistake.

The digital readout is easy to set. Simply rotate the plunger to the exact microliter setting desired. And you're ready to go.

Available from CMS, Fisher, A. H. Thomas, and other leading suppliers. Capacities from 10 to 100 μ l. Calibrated glass capillaries (50) and accessory maintenance kit included.

FOR MORE INFORMATION
CIRCLE THE NUMBER FOR THE NEW DRUMMOND PIPETTING EQUIPMENT CATALOG. TWELVE ILLUSTRATED PAGES THAT TELL ABOUT THE COMPLETE LINE OF DRUMMOND PIPETTING AIDS.

DRUMMOND SCIENTIFIC COMPANY

500 Parkway, Broomall, Pa. 19008

Call toll-free: (800) 523-7480
In Pa.: (215) 353-0200

For pipetting, Drummond took the work right out of your mouth.

